

# Post-operative oedema prevention

Preventing the build-up of post-operative oedema following orthopaedic surgery.



# A new approach

## Providing lower limb muscle contraction for the prevention and treatment of oedema

Easy-to-use, the geko™ is a battery powered, disposable neuromuscular electrostimulation device designed to increase blood flow in the deep veins of the leg.<sup>17</sup>

The geko™ device gently stimulates the common peroneal nerve contracting the calf and foot muscle pumps to prevent and treat oedema.

---

# 60%

The increase in blood flow is equal to 60%<sup>18</sup> of walking without a patient having to move.

---


# Zero

No wires or leads.  
Small, light and comfortable to wear.  
Silent in operation.

---

# 10g

Weighs just 10g.  
Quick and easy to fit.



# Swelling in the knee, in the acute recovery period following TKR\*, can impair quadriceps muscle strength and rehabilitation<sup>1,2,3</sup>

## Knee swelling after TKR can influence the post-operative results.<sup>2</sup>

- One month after TKR, quadriceps muscle strength drops 50% to 60% of pre-operative levels, despite the initiation of rehabilitation within 48 hours after surgery.<sup>4,5</sup>
- Muscle weakness, particularly in the quadriceps muscle, has profound functional consequences, especially in older individuals and can be associated with decreased gait speed, balance, stair-climbing ability, and ability to rise from a seated position, as well as an increased risk for falls.<sup>6,7</sup>
- Impaired functional performance can also be associated with length of stay and patient-reported outcomes<sup>8</sup> - and the associated knee swelling is known to increase rates of wound dehiscence and infection<sup>9</sup> and can delay rehabilitation.

## Causality:

Impairments in quadriceps muscle strength are caused by oedema.<sup>10</sup>

Spinal reflex activity from swelling can diminish quadriceps muscle strength. Over time, the fibres of the muscles can atrophy due to the lack of use.<sup>1</sup>

Strategies to address **early** quadriceps muscle weakness should address the underlying oedema.<sup>10</sup>

# The geko™ device is clinically proven to prevent the build-up of post-operative oedema<sup>11</sup>

- The geko™ device delivers mechanical compression by activating the calf and foot muscle pumps resulting in increased blood flow and the reduction of oedema.<sup>11</sup>
- The highly portable geko™ device can lower Ambulatory Venous Pressure (AVP) and Venous Transit Times (VTT) transferring tissue fluid back into the veins.<sup>12</sup>
- The geko™ device can prevent oedema and promote functional activity following foot surgery.<sup>13</sup>
- The geko™ device is well tolerated.<sup>14</sup>

## Associated geko™ benefits

The geko™ device also provides venous thrombosis prevention (VTE). NICE guidance (MTG19) supports use of the geko™ device for people who have a high risk of VTE and for whom pharmacological or other methods of VTE prevention are impractical or contraindicated.<sup>15,16</sup>




# Related geko™ evidence

An RCT comparing the effect of the geko™ device on post-operative oedema in Total Hip Replacement patients, demonstrates geko™ efficacy in preventing the build-up of post-operative oedema.<sup>13</sup>

The graph shows the change in knee circumference on the hip patient's operated leg.

## Change in knee circumference, operated leg



# Providing lower limb muscle contraction for the prevention and treatment of oedema

## Self-contained and wearable, the geko™ device is:

- Simple and easy to use.
- Small and light (weighing just 10g) with no wires or leads, enables the patient to be as mobile as possible.

### References

1. Stevens-Lapsley E et al. Early Neuromuscular Electrical Stimulation to Improve Quadriceps Muscle Strength After Total Knee Arthroplasty: A Randomized Controlled Trial. *Phys Ther* 2012 Feb; 92(2): 210–226. Published online 2011 Nov 17. doi: [10.2522/ptj.20110124].
2. Pua YH. The Time Course of Knee Swelling Post Total Knee Arthroplasty and Its Associations with Quadriceps Strength and Gait Speed. *J Arthroplasty*. 2015 Jul;30(7):1215–9. doi: 10.1016/j.arth.2015.02.010. Epub 2015 Feb 19.
3. Palmieri-Smith R M et al. Pain and Effusion and Quadriceps Activation and Strength. *J Athl Train*. 2013 Mar-Apr; 48(2): 186–191. doi: [10.4085/1062-6050-48.2.10].
4. Milzner RL et al. Early Quadriceps Strength Loss After Total Knee Arthroplasty. The Contributions of Muscle Atrophy and Failure of Voluntary Muscle Activation. *J Bone Joint Surg Am*. 2005 May;87(5):1047–53. <https://www.ncbi.nlm.nih.gov/pubmed/15866968>
5. Bade MJ et al. Outcomes Before and After Total Knee Arthroplasty Compared to Healthy Adults. *J Orthop Sports Phys Ther*. 2010 Sep;40(9):559–67. doi: 10.2519/jospt.2010.3317.
6. Rantanen T et al. Association of Muscle Strength With Maximum Walking Speed in Disabled Older Women. *Am J Phys Med Rehabil*. 1998 Jul-Aug;77(4):299–305.
7. Milzner RL et al. Quadriceps Strength and the Time Course of Functional Recovery After Total Knee Arthroplasty. *J Orthop Sports Phys Ther*. 2005 Jul;35(7):424–36.
8. Brock TM et al. Short-stretch Inelastic Compression Bandage in Knee Swelling Following Total Knee Arthroplasty Study (STICKS): Study Protocol for a Randomised Controlled Feasibility Study. *Trials* 2015;16:87.
9. Yu GV et al. The Jones Compression Bandage. Review and Clinical Applications. *J Am Podiatr Med Assoc* 2002;92:221–31.
10. Holm B et al. Loss of Knee-Extension Strength is Related to Knee Swelling after Total Knee Arthroplasty. *Arch Phys Med Rehabil*. 2010 Nov;91(11):1770–6. doi: 10.1016/j.apmr.2010.07.229.
11. Wainwright TW et al. A Feasibility Randomised Controlled Trial to Evaluate the Effectiveness of a Novel Neuromuscular Electro-stimulation Device in Preventing the Formation of Oedema Following Total Hip Replacement Surgery. *Heliyon* 18 Jul 2018- Volume 4, Issue 7.
12. M Khanbhai, Vhansrani, J Sultan, J Burke, CN McCollum. The Effect of Neuromuscular Electrostimulation on Lower Limb Venous Physiology. Academic Surgery Unit, Institute of Cardiovascular Sciences, Manchester Academic Health Science Centre. SARS Annual meeting, 2014. University of Cambridge, 8th & 9th January 2014.
13. Case studies. Data on file Firstkind Limited.
14. P Baker et al. Abstract 2017. Data on file Firstkind Limited: The geko™ Neuromuscular Electrostimulation Device Reduces Pre-Operative Oedema and Accelerates Readiness to Theatre in Patients Requiring Open Reduction Internal Fixation for Ankle Fracture.
15. NICE Guidance (MTG19). Published June 2014.
16. Yilmaz S et al. Potential Role of Electrostimulation in Augmentation of Venous Blood Flow After Total Knee Replacement: A Pilot Study. *Phlebology*. 2016 May;31(4):251–6. doi: 10.1177/0268355515580473. Epub 2015 Apr 6.
17. Nicolaidas A & Griffin M. Measurement of Blood Flow in the Deep Veins of the Lower Limb Using the geko™ Neuromuscular Electro-stimulation Device. *International Angiology* 2016 August;35(4):406–10.
18. Tucker A et al. Augmentation of Venous, Arterial and Microvascular Blood Supply in the Leg by Isometric Neuromuscular Stimulation via the Peroneal Nerve. *Int J Angiol*. 2010 Spring; 19(1): e31–e37.
19. Klabunde, RE (2014). Cardiovascular Physiology Concepts. Available at <http://cvphysiology.com/Microcirculation/M010> [Accessed 21 Feb. 2018].



**firstkind**  
living science

**Firstkind Ltd, Hawk House, Peregrine Business Park, High Wycombe, Buckinghamshire, HP13 7DL, United Kingdom.**

**T: +44 (0)845 2222 921  
W: [www.gekodevices.com](http://www.gekodevices.com)**

MPALLORTH0401